



# State of New Hampshire

## Department of Safety



Division of Fire Standards and Training & Emergency Medical Services  
 Richard M. Flynn Fire Academy  
 98 Smokey Bear Boulevard, Concord, New Hampshire  
 Mailing Address: 33 Hazen Drive, Concord, New Hampshire 03305-0002

John J. Barthelmes  
 Commissioner

Robert Quinn and Richard C. Bailey, Jr.  
 Assistant Commissioners

Deborah A. Pendergast  
 Director

# CLINICAL BULLETIN

Bulletin #	Title			Date Issued
49	<b>LIFEPAK® 15 Monitor/Defibrillator by Stryker: Recall</b>			February 4, 2019
Superseded	Released By	Approved By	Source	
	Vicki Blanchard, NRP	Joey Scollan, MD	FDA	

A new MedWatch Safety Alert recall: **LIFEPAK® 15 Monitor/Defibrillator by Stryker: Recall - Due to Lock Up (blank monitor display with LED lights on) after a Defibrillation Shock is Delivered** was just added to the FDA Recalls web page.

You can read the complete recall at: <https://www.fda.gov/Safety/Recalls/ucm630455.htm>